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1. Akino Y, Tohyama N, Akita K, et al. Modalities and techniques used for stereotactic radiotherapy, intensity-modulated radiotherapy, and image-guided radiotherapy: A 2018 survey by the Japan Society of Medical Physics. *Phys Med* 2009; 64: 182-187
2. 放射線治療品質管理機構：放射線治療装置導入に関するコミッショニング必要期間について。2008。  
[https://www.qcrt.org/common/pdf/commissioning\\_kikan.pdf](https://www.qcrt.org/common/pdf/commissioning_kikan.pdf), (2 Sep 2022, date last accessed).
3. Klein EE, Hanley J, Bayouth J, et al. Task Group 142 report: quality assurance of medical accelerators. *Med Phys* 2009; 36: 4197-4212.
4. Das IJ, Cheng CW, Watts RJ, et al. Accelerator beam data commissioning equipment and procedures: report of the TG-106 of the Therapy Physics Committee of the AAPM. *Med Phys* 2008; 35: 4186-4215.
5. Isono M, Akino Y, Mizuno H, et al. Inter-unit variability of multi-leaf collimator parameters for IMRT and VMAT treatment planning: a multi-institutional survey. *J Radiat Res* 2020; 61: 307-313.
6. Smilowitz JB, Das IJ, Feygelman V, et al. AAPM Medical Physics Practice Guideline 5.b.: Commissioning and QA of Treatment Planning Dose Calculations - Megavoltage Photon and Electron Beams. *J Appl Clin Med Phys* 2022; e13641.
7. Miften M, Olch A, Mihailidis D, et al. Tolerance limits and methodologies for IMRT measurement-based verification QA: Recommendations of AAPM Task Group No. 218. *Med Phys* 2018; 45: e53-e83.
8. International Atomic Energy Agency, Commissioning, and quality assurance of computerized planning systems for radiation treatment of cancer, Technical Report Series No. 430, IAEA, Vienna, 2004
9. Huq MS, Fraass BA, Dunscombe PB, et al. The report of Task Group 100 of the AAPM: Application of risk analysis methods to radiation therapy quality management. *Med Phys* 2016; 43: 4209-4261.
10. 安全は偶然に生まれない 放射線腫瘍学の品質とケアを向上させるための枠組み 日本語訳  
[https://www.jastro.or.jp/medicalpersonnel/safety/safety\\_is\\_no\\_accident.pdf](https://www.jastro.or.jp/medicalpersonnel/safety/safety_is_no_accident.pdf), (2 Sep 2022, date last accessed).

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1. Klein EE, Hanley J, Bayouth J, et al. Task Group 142 report: Quality assurance of medical accelerators. *Med Phys* 2009; 36: 4197-4212.
2. Hanley J, Dresser S, Simon W, et al. AAPM Task Group 198 Report: An implementation guide for TG 142 quality assurance of medical

accelerators. *Med Phys* 2021; 48: e830-e885.

3. 日本放射線腫瘍学会 QA 委員会. 外部放射線治療における Quality Assurance (QA) システムガイドライン 2016 年版, 金原出版, 東京, 2016.
4. 黒岡将彦. Linac・RTPS の QA/QC (コミッションング). 日放線技会放射線治療分科会誌 2009; 23: 8-18.
5. Langen KM, Papanikolaou N, Balog J, et al. QA for helical tomotherapy: report of the AAPM Task Group 148. *Med Phys* 2021; 37: 4817-4853.
6. Dieterich S, Cavedon C, Chuang CF, et al. Report of AAPM TG 135: quality assurance for robotic radiosurgery. *Med Phys* 2011; 38: 2914-2936.
7. Xiao Y, Kry SF, Popple R, et al. Flattening filter-free accelerators: a report from the AAPM Therapy Emerging Technology Assessment Work Group. *J Appl Clin Med Phys* 2015; 16: 12-29.
8. Otto K. Volumetric modulated arc therapy: IMRT in a single gantry arc. *Med Phys* 2008; 35: 310-317.
9. Barish RJ, Fleischman RC and Pipman YM. Teletherapy beam characteristics: The first second. *Med Phys* 1987; 14:657-661.
10. Sonke JJ, Brand B and van Herk M. Focal spot motion of linear accelerators and its effect on portal image analysis. *Med Phys* 2003; 30: 1067-1075.
11. Waldron TJ: Functional requirements for IMRT. In: Palta JR and Mackie TR eds. Intensity-Modulated Radiation Therapy: The State of the Art. Wisconsin, Medical Physics Publishing, 2003, 373-400.
12. Ezzell GA and Chungbin S. The overshoot phenomenon in step-and-shoot IMRT delivery. *J Appl Clin Med Phys* 2: 138-148, 2001.
13. Li J, Wiersma RD, Stepaniak CJ, et al: Improvements in dose accuracy delivered with static-MLC IMRT on an integrated linear accelerator control system. *Med Phys* 2012; 39: 2456-2462.
14. Yu CY, Wan SW, Weng YC, et al. Three corrections for overshoot effect improved the dose for step-and-shoot intensity-modulated radiation therapy. *PLoS One* 2021; 16: e0250243.
15. Litzenberg DW, Moran JM and Fraass BA. Incorporation of realistic delivery limitations into dynamic MLC treatment delivery. *Med Phys* 2002; 29: 810-820.
16. Xia P, Chuang CF and Verhey LJ. Communication and sampling rate limitations in IMRT delivery with a dynamic multileaf collimator system. *Med Phys* 2002; 29: 412-423.
17. Low DA, Sohn JW, Klein EE, et al. Characterization of a commercial multileaf collimator used for intensity modulated radiation therapy. *Med Phys* 2002; 28: 752-756.
18. Kerns JR, Childress N, Kry SF. A multi-institution evaluation of MLC log files and performance in IMRT delivery. *Radiat Oncol* 9:176,2014.
19. LoSasso T, Chui CS and Ling C: Physical and dosimetric aspects of a multileaf collimation system used in the dynamic mode for implementing intensity modulated radiotherapy. *Med Phys* 1998; 25: 1919-1927.
20. Arnfield MR, Siebers JV, Kim JO, et al. A method for determining multileaf collimator transmission and scatter for dynamic intensity modulated radiotherapy. *Med Phys* 2000; 27: 2231-2241.
21. Lim TY, Dragojević I, Hoffman D, et al. Characterization of the Halcyon™ multileaf collimator system. *J Appl Clin Med Phys* 2019; 20: 106-114.
22. Agnew A, Agnew CE, Grattan MW, et al. Monitoring daily MLC positional errors using trajectory log files and EPID measurements for IMRT and VMAT deliveries. *Phys Med Biol* 2014; 59: N49-63.
23. Palta JR, Kim S, Jonathan GL, et al. Tolerance Limits and action levels for planning and delivery of IMRT. In: Palta JR and Mackie TR eds. Intensity-Modulated Radiation Therapy: The State of the Art. Wisconsin, Medical Physics Publishing, 2003, 593-612.
24. Mu G, Ludlum E, Xia P. Impact of MLC leaf position errors on simple and complex IMRT plans for head and neck cancer. *Phys Med Biol* 2008; 53: 77-88.
25. Rangel A, Dunscombe P. Tolerances on MLC leaf position accuracy for IMRT delivery with a dynamic MLC. *Med Phys* 2009; 36: 3304-3309.
26. Oliver M, Gagne I, Bush K, et al. Clinical significance of multi-leaf collimator positional errors for volumetric modulated arc therapy. *Radiother*

*Oncol* 2010; 97: 554-560.

27. Tatsumi D, Hosono MN, Nakada R, et al. Direct impact analysis of multi-leaf collimator leaf position errors on dose distributions in volumetric modulated arc therapy: a pass rate calculation between measured planar doses with and without the position errors. *Phys Med Biol* 2011; 56: N237-46.
28. LoSasso T, Chui CS and Ling CC. Comprehensive quality assurance for the delivery of intensity modulated radiotherapy with a multileaf collimator used in the dynamic mode. *Med Phys* 2001; 28: 2209-2219.
29. Wen N, Zhao B, Kim J, et al. IMRT and RapidArc commissioning of a TrueBeam linear accelerator using TG-119 protocol cases. *J Appl Clin Med Phys* 2014; 15: 74-88.
30. LoSasso T. Acceptance testing and commissioning of IMRT. In: Fuks Z, Leibel SA, Ling CC eds. A practical guide to intensity-modulated radiation therapy. Wisconsin, Medical Physics Publishing, 2003, 123-146.
31. Bedford JL and Warrington A. Commissioning of volumetric modulated arc therapy (VMAT). *Int J Radiat Oncol Biol Phys* 2009; 73: 537-545.
32. Ling CC, Zhang P, Archambault Y, et al. Commissioning and quality assurance of RapidArc radiotherapy delivery system. *Int J Radiat Oncol Biol Phys* 2008; 72: 575-581.
33. Haga A, Nakagawa K, Shiraiishi K, et al. Quality assurance of volumetric modulated arc therapy using Elekta Synergy. *Acta Oncol* 2009; 48: 1193-1197.
34. Schreibmann E, Dhabaan A, Elder E, et al. Patient-specific quality assurance method for VMAT treatment delivery. *Med Phys* 2009; 36: 4530-4535.
35. トモセラピー物理ワーキンググループ. トモセラピーシステムの物理・技術的 Q&A. 日立メディコ株式会社, 2014.
36. Accuray. Radixact® Treatment Delivery System, Physics Essentials Guide. 1067608-ENG A, 2021.
37. Fenwick JD, Tome WA, Jaradat HA, et al. Quality assurance of a helical tomotherapy machine. *Phys Med Biol* 2004; 49: 2933-2953.
38. Balog J and Mackie TR. Benchmarking beam alignment for a clinical helical tomotherapy device. *Med Phys* 2003; 30: 1118-1127.
39. Balog J, Olivera G, and Kapatoes J. Clinical helical tomotherapy commissioning dosimetry. *Med Phys* 2003; 30: 3097-3106.
40. 日本医学物理学会編. 外部放射線治療における水吸収線量の標準計測法 (標準計測法 12) . 通商産業研究社, 東京, 2012.
41. Schnarr E, Beneke M, Casey D, et al. Feasibility of real-time motion management with helical tomotherapy. *Med Phys* 2018; 45: 1329-1337.
42. Chen GP, Tai A, Keiper TD, et al. Technical Note: Comprehensive performance tests of the first clinical real-time motion tracking and compensation system using MLC and jaws. *Med Phys* 2020; 47: 2814-2825.
43. Ferris WS, Kissick MW, Bayouth JE, et al. Evaluation of radixact motion synchrony for 3D respiratory motion: Modeling accuracy and dosimetric fidelity. *J Appl Clin Med Phys* 2020; 21:96-106.
44. Goddard L, Jeong K, Tomé WA. Commissioning and routine quality assurance of the Radixact Synchrony system. *Med Phys* 2022; 49: 1181-1195.
45. Acuray. CyberKnife®治療実行システム, フィジックス基本ガイド, 1062939-JPN A, 2018.
46. サイバーナイフ QA 委員会. サイバーナイフ QA マニュアルII, 2018.
47. LoSasso T. IMRT delivery system QA. In: Palta J and Mackie TR, editors. Intensity modulated radiation therapy: the state of the art. Madison, WI: Medical Physics Publishing; 2003, 561-591.
48. Bayouth JE, Wendt D, Morrill SM. MLC quality assurance techniques for IMRT applications. *Med Phys* 2003; 30: 743-750.

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26. Das IJ, Ding GX and Ahnesjo A. Small fields: nonequilibrium radiation dosimetry. *Med Phys* 2008; 35: 206-215.
27. Kawachi T, Saitoh H, Inoue M, et al. Reference dosimetry condition and beam quality correction factor for CyberKnife beam. *Med Phys* 2008; 35: 4591-4598.
28. Jordan TJ and Williams PC. The design and performance characteristics of a multileaf collimator. *Phys Med Biol* 1994; 39: 231-251.
29. Yan G, Fox C, Liu C, et al. The extraction of true profiles for TPS commissioning and its impact on IMRT patient-specific QA. *Med Phys* 2008; 35: 3661-3670.
30. LoSasso T, Chui CS, and Ling CC. Physical and dosimetric aspects of a multileaf collimation system used in the dynamic mode for implementing intensity modulated radiotherapy. *Med Phys* 1998; 25: 1919-1927.
31. Arnfield MR, Siebers JV, Kim JO, et al. A method for determining multileaf collimator transmission and scatter for dynamic intensity modulated radiotherapy. *Med Phys* 2000; 27: 2231-2241.
32. 岸和馬. DMLC による IMRT : IMRT における QA と線量照合. 日放線技会放射線治療分科会誌 2002; 16: 44-50.
33. Isono M, Akino Y, Mizuno H, et al. Inter-unit variability of multi-leaf collimator parameters for IMRT and VMAT treatment planning: a multi-institutional survey. *J Radiat Res* 2020; 61: 307-313.
34. Kim J, Han JS, Hsia AT, et al. Relationship between dosimetric leaf gap and dose calculation errors for high definition multi-leaf collimators in radiotherapy. *Phys Imaging Radiat Oncol.* 2018; 5: 31-36.
35. Kielar KN, Mok E, Hsu A, et al. Verification of dosimetric accuracy on the TrueBeam STx: rounded leaf effect of the high definition MLC, *Med Phys.* 2012; 39: 6360-71.
36. Szpala S, Cao F, Kohli K. On using the dosimetric leaf gap to model the rounded leaf ends in VMAT/RapidArc plans. *J Appl Clin Med Phys* 2014; 15: 67-84.
37. Wen N, Zhao B, Kim J, et al. IMRT and RapidArc commissioning of a TrueBeam linear accelerator using TG-119 protocol cases. *J Appl Clin Med Phys* 2014; 15: 74-88.
38. Tsuruta Y, Nakata M, Nakamura M, et al. Dosimetric comparison of Acuros XB, AAA, and XVMC in stereotactic body radiotherapy for lung cancer. *Med Phys* 2014; 41: 081715.
39. Kry SF, Feygelman V, Balter P, et al. AAPM Task Group 329: Reference dose specification for dose calculations: Dose-to-water or dose-to-muscle? *Med Phys* 2020; 47: e52-e64.
40. Anetai Y, Sumida I, Kumazaki Y, et al. Assessment of using a gamma index analysis for patient-specific quality assurance in Japan. *J Appl Clin Med Phys* 2022; 26: e13745.

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1. Hanley J, Dresser S, Simon W, et al. AAPM Task Group 198 Report: An implementation guide for TG 142 quality assurance of medical accelerators. *Med Phys* 2021; 48: e830-e885.
2. Mutic S, Palta JR, Butker EK, et al. Quality assurance for computed-tomography simulators and the computed-tomography-simulation process: report of the AAPM Radiation Therapy Committee Task Group No. 66. *Med Phys* 2003; 30: 2762-2792.
3. Li H, Noel C, Chen H, et al. Clinical evaluation of a commercial orthopedic metal artifact reduction tool for CT simulations in radiation therapy. *Med Phys* 2012; 39: 7507-7517.
4. White G, Wilson I. Photon, electron, proton and neutron interaction data for body tissues. *ICRU report* 1992; 46: 7-200.
5. INTERNATIONAL ATOMIC ENERGY AGENCY. Commissioning and Quality Assurance of Computerized Planning Systems for Radiation Treatment of Cancer. Vienna, 2004.
6. Smilowitz JB, Das IJ, Feygelman V, et al. AAPM Medical Physics Practice Guideline 5.a.: Commissioning and QA of Treatment Planning Dose Calculations — Megavoltage Photon and Electron Beams. *J Appl Clin Med Phys* 2015; 16: 14-34.

7. Davis AT, Palmer AL, Nisbet A. Can CT scan protocols used for radiotherapy treatment planning be adjusted to optimize image quality and patient dose? A systematic review. *Br J Radiol* 2017; 90: 20160406.
8. Sadagopan R, Bencomo JA, Martin RL, et al. Characterization and clinical evaluation of a novel IMRT quality assurance system. *J Appl Clin Med Phys* 2009; 10: 104-119.
9. Tyagi N, Fontenla S, Zhang J, et al. Dosimetric and workflow evaluation of first commercial synthetic CT software for clinical use in pelvis. *Phys Med Biol*. 2017; 62: 2961-2975.
10. 日本医学物理学会, 日本高精度放射線外部照射研究会, 日本放射線技術学会, 他. 呼吸性移動対策を伴う放射線治療に関するガイドライン 2019.
11. Jiang SB, Wolfgang J, Mageras GS. Quality assurance challenges for motion-adaptive radiation therapy: gating, breath holding, and four-dimensional computed tomography. *Int J Radiat Oncol Biol Phys* 2008; 71: S103-S107.
12. Lambrecht M, Sonke JJ, Nestle U, et al. Quality assurance of four-dimensional computed tomography in a multicentre trial of stereotactic body radiotherapy of centrally located lung tumours. *Phys Imaging Radiat Oncol* 2018; 8: 57-62.
13. Yamamoto T, Langner U, Loo BW Jr., et al. Retrospective analysis of artifacts in four-dimensional CT images of 50 abdominal and thoracic radiotherapy patients. *Int J Radiat Oncol Biol Phys* 2008; 72: 1250-1258.
14. Nakamura M, Narita Y, Sawada A, et al. Impact of motion velocity on four-dimensional target volumes: a phantom study. *Med Phys* 2009; 36: 1610-1617.
15. 厚生労働省. 医療法施行規則の一部を改正する省令の施行等について, 医政発 0312 第 7 号, 平成 31 年 3 月 12 日.
16. Ellingson BM, Bendszus M, Boxerman J, et al. Consensus recommendations for a standardized Brain Tumor Imaging Protocol in clinical trials. *Neuro Oncol* 2015; 17: 1188-1198.
17. INTERNATIONAL ATOMIC ENERGY AGENCY. The role of PET/CT in radiation treatment planning for cancer patient treatment. Vienna, 2008.
18. Paulson ES, Erickson B, Schultz C, et al. Comprehensive MRI simulation methodology using a dedicated MRI scanner in radiation oncology for external beam radiation treatment planning. *Med Phys* 2015; 42: 28-39.
19. Glide-Hurst CK, Paulson ES, McGee, K, et al. Task group 284 report: magnetic resonance imaging simulation in radiotherapy: considerations for clinical implementation, optimization, and quality assurance. *Med Phys* 2021; 48: e636-e670.
20. Konert T, Vogel W, MacManus MP, et al. PET/CT imaging for target volume delineation in curative intent radiotherapy of non-small cell lung cancer: IAEA consensus report 2014. *Radiother Oncol* 2015; 116: 27-34.
21. Kavaluus H, Nousiainen K, Kaijaluoto S, et al. Determination of acceptance criteria for geometric accuracy of magnetic resonance imaging scanners used in radiotherapy planning. *Phys Imaging Radiat Oncol* 2021; 17: 58-64.
22. 日本放射線腫瘍学会. 放射線治療における非剛体画像レジストレーション利用のためのガイドライン 2018 年版.
23. 日本放射線腫瘍学会. 画像誘導放射線治療の臨床施行のためのガイドライン 2022.
24. Posiewnik M, Piotrowski T. A review of cone-beam CT applications for adaptive radiotherapy of prostate cancer. *Phys Med*. 2019; 59: 13-21.
25. Guan H, Dong H. Dose calculation accuracy using cone-beam CT (CBCT) for pelvic adaptive radiotherapy. *Phys Med Biol*. 2009; 54(20): 6239-6250.
26. Bissonnette JP, Balter PA, Dong L, et al. Quality assurance for image-guided radiation therapy utilizing CT-based technologies: a report of the AAPM TG-179. *Med Phys*. 2012; 39(4): 1946-1963.
27. Jarema T, Aland T. Using the iterative kV CBCT reconstruction on the Varian Halcyon linear accelerator for radiation therapy planning for pelvis patients. *Phys Med*. 2019; 68: 112-116.
28. Langen KM, Papanikolaou N, Balog J, et al. QA for helical tomotherapy: report of the AAPM Task Group 148. *Med Phys* 2010; 37: 4817-4853.

29. Petit SF, van Elmpt WJ, Lambin P, et al. Dose recalculation in megavoltage cone-beam CT for treatment evaluation: removal of cupping and truncation artefacts in scans of the thorax and abdomen. *Radiother Oncol* 2010; 94: 359-366.
30. Ding GX, Alaei P, Curran B, et al. Image guidance doses delivered during radiotherapy: Quantification, management, and reduction: Report of the AAPM Therapy Physics Committee Task Group 180. *Med Phys* 2018; 45: e84-e99.
31. Kutcher GJ, Coia L, Gillin M, et al. Comprehensive QA for radiation oncology: report of AAPM Radiation Therapy Committee Task Group 40. *Med Phys* 1994; 21: 581-618.
32. Hynds S, McGarry CK, Mitchell DM, et al. Assessing the daily consistency of bladder filling using an ultrasonic Bladderscan device in men receiving radical conformal radiotherapy for prostate cancer. *Br J Radiol.* 2011; 84: 813-818.
33. Shirato H, Suzuki K, Sharp GC, et al. Speed and amplitude of lung tumor motion precisely detected in four-dimensional setup and in real-time tumor-tracking radiotherapy. *Int J Radiat Oncol Biol Phys* 2006; 64: 1229-1236.
34. Shiinoki T, Kawamura S, Uehara T, et al. Evaluation of a combined respiratory-gating system comprising the TrueBeam linear accelerator and a new real-time tumor-tracking radiotherapy system: a preliminary study. *J Appl Clin Med Phys* 2016; 17: 202-213.
35. Hoogeman M, Prevost JB, Nuytens J, et al. Clinical accuracy of the respiratory tumor tracking system of the cyberknife: assessment by analysis of log files. *Int J Radiat Oncol Biol Phys* 2009; 74: 297-303.
36. Mukumoto N, Nakamura M, Yamada M, et al. Intrafractional tracking accuracy in infrared marker-based hybrid dynamic tumour-tracking irradiation with a gimballed linac. *Radiother Oncol* 2014; 111: 301-305.
37. Akimoto M, Nakamura M, Mukumoto N, et al. Predictive uncertainty in infrared marker-based dynamic tumor tracking with Vero4DRT. *Med Phys* 2013; 40: 091705.

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1. Srivastava SP, Cheng CW, Das IJ. The effect of slice thickness on target and organs at risk volumes, dosimetric coverage and radiobiological impact in IMRT planning. *Clin Transl Oncol* 2016; 18: 469-79.
2. Anetai Y, Sumida I, Kumazaki Y, et al. Assessment of using a gamma index analysis for patient-specific quality assurance in Japan. *J Appl Clin Med Phys* 2022; 26: e13745.
3. 日本放射線腫瘍学会他. 放射線治療における非剛体画像レジストレーション利用のためのガイドライン 2018年版.
4. International Commission on Radiation Units and Measurements: Prescribing, Recording and Reporting photon Beam Therapy, ICRU Report 50. Bethesda, ICRU Publications, 1993.
5. International Commission on Radiation Units and Measurements: Prescribing, Recording and Reporting Photon Beam Therapy (Supplement to ICRU Report 50), ICRU Report 62. Bethesda, ICRU Publications, 1999.
6. International Commission on Radiation Units and Measurements: Prescribing, Recording and Reporting Photon-Beam Intensity-Modulated Radiation Therapy (IMRT), ICRU Report 83. Bethesda, ICRU Publications, 2010.
7. Mayo CS, Moran JM, Bosch W, et al. American Association of Physicists in Medicine Task Group 263: Standardizing Nomenclatures in Radiation Oncology. *Int J Radiat Oncol Biol Phys.* 2018; 100: 1057-1066.
8. van Herk M et al. The probability of correct target dosage: dose-population histograms for deriving treatment margins in radiotherapy. *Int J Radiat Oncol Biol Phys* 2000; 47: 1121-1135.
9. Olch AJ, Gerig L, Li H, et al. Dosimetric effects caused by couch tops and immobilization devices: report of AAPM Task Group 176. *Med Phys* 2014; 41: 061501.
10. Giorgia N, Antonella F, Alessandro C, et al. Planning strategies in volumetric modulated arc therapy for breast. *Med Phys* 2011; 38: 4025-4031.
11. Tyran M, Tallet A, Resbeut M, et al. Safety and benefit of using a virtual bolus during treatment planning for breast cancer treated with arc therapy. *J Appl Clin Med Phys* 2018; 19: 463-472.

12. Ong CL, Dahele M, Slotman BJ, et al. Dosimetric impact of the interplay effect during stereotactic lung radiation therapy delivery using flattening filter-free beams and volumetric modulated arc therapy. *Int J Radiat Oncol Biol Phys* 2013; 86: 743-748.
13. Court LE, Seco J, Lu XQ, et al. Use of a realistic breathing lung phantom to evaluate dose delivery errors. *Med Phys* 2010; 37: 5850-5857.
14. Court L, Wagar M, Bogdanov M, et al. Use of reduced dose rate when treating moving tumors using dynamic IMRT. *J Appl Clin Med Phys* 2010; 12: 3276.
15. Netherton T, Li Y, Nitsch P, et al. Interplay effect on a 6-MV flattening-filter-free linear accelerator with high dose rate and fast multi-leaf collimator motion treating breast and lung phantoms. *Med Phys* 2018; 45: 2369-2376.
16. Huq MS, Fraass BA, Dunscombe PB, et al. The report of Task Group 100 of the AAPM: Application of risk analysis methods to radiation therapy quality management. *Med Phys* 2016; 43: 4209.
17. Ford E, Conroy L, Dong L, et al. Strategies for effective physics plan and chart review in radiation therapy: Report of AAPM Task Group 275. *Med Phys* 2020; 47: e236-e272.
18. Xia P, Sintay BJ, Colussi VC, et al. Medical Physics Practice Guideline (MPPG) 11.a: Plan and chart review in external beam radiotherapy and brachytherapy. *J Appl Clin Med Phys* 2021; 22: 4-19.
19. 日本放射線腫瘍学会. 外部放射線治療における QA システムガイドライン 2016 年版, 金原出版, 東京, 2016.
20. 厚生労働省. 医療情報システムの安全管理に関するガイドライン第 5.2 版. 令和 4 年 3 月.
21. 日本放射線技術学会. 画像情報の確定に関するガイドライン Ver.2.1. 平成 26 年 8 月 18 日.
22. Siochi RA, Balter P, Bloch CD, et al. Report of Task Group 201 of the American Association of Physicists in Medicine: Quality management of external beam therapy data transfer. *Med Phys* 2021; 48: e86-e114.
23. Vrtovec T, Močnik D, Strojani P, et al. Auto-segmentation of organs at risk for head and neck radiotherapy planning: From atlas-based to deep learning methods. *Med Phys* 2020; 47: e929-e950.
24. Kraysenbuehl J, Zamburlini M, Ghandour S, et al. Planning comparison of five automated treatment planning solutions for locally advanced head and neck cancer [published correction appears in *Radiat Oncol* 2018; 13: 225.]. *Radiat Oncol*. 2018; 13: 170.
25. Babier A, Mahmood R, McNiven AL, et al. Knowledge-based automated planning with three-dimensional generative adversarial networks. *Med Phys* 2020; 47: 297-306.
26. Sherer MV, Lin D, Elguindi S, et al. Metrics to evaluate the performance of auto-segmentation for radiation treatment planning: A critical review. *Radiother Oncol* 2021; 160: 185-191.
27. Chang ATY, Hung AWM, Cheung FWK, et al. Comparison of Planning Quality and Efficiency Between Conventional and Knowledge-based Algorithms in Nasopharyngeal Cancer Patients Using Intensity Modulated Radiation Therapy. *Int J Radiat Oncol Biol Phys* 2016; 95: 981-990.
28. Kubo K, Monzen H, Ishii K, et al. Dosimetric comparison of RapidPlan and manually optimized plans in volumetric modulated arc therapy for prostate cancer. *Phys Med* 2017; 44: 199-204.
29. Wall PDH, Fontenot JD. Evaluation of complexity and deliverability of prostate cancer treatment plans designed with a knowledge-based VMAT planning technique. *J Appl Clin Med Phys* 2020; 21: 69-77.

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1. Ezzell GA, Galvin JM, Low D, et al. Guidance document on delivery, treatment planning, and clinical implementation of IMRT: report of the IMRT subcommittee of the AAPM radiation therapy committee. *Med Phys* 2003; 30: 2089-2115.
2. Ezzell GA, Burmeister JW, Dogan N, et al. IMRT commissioning: multiple institution planning and dosimetry comparisons, a report from AAPM Task Group 119. *Med Phys* 2000; 36: 5359-5373.
3. Moran JM, Dempsey M, Eisbruch A, et al. Safety considerations for IMRT: executive summary. *Med Phys* 2011; 38: 5067.
4. Hartford AC, Galvin JM, Beyer DC, et al. American College of Radiology (ACR) and American Society for Radiation Oncology (ASTRO)

practice guideline for intensity-modulated radiation therapy (IMRT). *Am J Clin Oncol* 2012; 35: 612–617.

5. Bogdanich W. Radiation Offers New Cures, and Ways to Do Harm. New York: The New York Times, 2010.
6. Bogdanich W. As Technology Surges, Radiation Safeguards Lag. New York: The New York Times, 2010.
7. Miften M, Olch A, Mihailidis D, et al. Tolerance limits and methodologies for IMRT measurement-based verification QA: Recommendations of AAPM Task Group No. 218. *Med Phys* 2018; 45: e53–e83.
8. Zhu TC, Stathakis S, Clark JR, et al: Report of AAPM Task Group 219 on independent calculation-based dose/MU verification for IMRT. *Med Phys* 2021; 48: e808–e829.
9. McEwen M, DeWerd L, Ibbott G, et al. Addendum to the AAPM's TG-51 protocol for clinical reference dosimetry of high energy photon beams. *Med Phys* 2014; 41: 041501.
10. Palmans H, Andreo P, Huq MS, et al. Dosimetry of small static fields used in external beam radiotherapy: an international code of practice for reference and relative dose determination. IAEA Technical Report Series No. 483. IAEA, Vienna, 2017.
11. Low DA, Moran JM, Dempsey JF, et al. Dosimetry tools and techniques for IMRT. *Med Phys* 2011; 38: 1313–1338.
12. Anetai Y, Sumida I, Kumazaki Y, et al. Assessment of using a gamma index analysis for patient-specific quality assurance in Japan. *J Appl Clin Med Phys* 2022; 26: e13745.
13. Suchowerska N, Hoban P, Butson M, et al. Directional dependence in film dosimetry: radiographic film and radiochromic film. *Phys Med Biol* 2001; 46: 1391–1397.
14. Pai S, Das IJ, Dempsey JF, et al. TG-69: Radiographic film for megavoltage beam dosimetry. *Med Phys* 2007; 34: 2228–2258.
15. Greer PB and Popescu CC. Dosimetric properties of an amorphous silicon electronic portal imaging device for verification of dynamic intensity modulated radiation therapy. *Med Phys* 2003; 30: 1618–1627.
16. Anjum MN, Parker W, Ruo R, et al. IMRT quality assurance using a second treatment planning system. *Med Dosim* 2010; 35: 274–279.
17. Fontenot JD. Evaluation of a novel secondary check tool for intensity-modulated radiotherapy treatment planning. *J Appl Clin Med Phys* 2014; 15: 4990.
18. Tsuruta Y, Nakamura M, Miyabe Y, et al. Use of a second-dose calculation algorithm to check dosimetric parameters for the dose distribution of a first-dose calculation algorithm for lung SBRT plans. *Phys Med* 2017; 44: 86–95.
19. Ono T, Mitsuyoshi T, Shintani T, et al: Independent calculation-based verification of volumetric-modulated arc therapy-stereotactic body radiotherapy plans for lung cancer. *J Appl Clin Med Phys* 2020; 21: 135–143.
20. McDonald DG, Jacqmin DJ, Mart CJ, et al. Validation of a modern second-check dosimetry system using a novel verification phantom. *J Appl Clin Med Phys* 2017; 18: 170–177.
21. Bresciani S, Poli M, Miranti A, et al. Comparison of two different EPID-based solutions performing pretreatment quality assurance: 2D portal dosimetry versus 3D forward projection method. *Phys Med* 2018; 52: 65–71.
22. Van Dyk J, Barnett RB, Cygler JE, et al. Commissioning and quality assurance of treatment planning computers. *Int J Radiat Oncol Biol Phys* 1993; 26: 261–273.
23. Harms WB Sr, Low DA, Wong JW, et al. A software tool for the quantitative evaluation of 3D dose calculation algorithms. *Med Phys* 1998; 25: 1830–1836.
24. Low DA, Harms WB, Mutic S, et al. A technique for the quantitative evaluation of dose distributions. *Med Phys* 1998; 25: 656–661.
25. Low DA, Dempsey JF. Evaluation of the gamma dose distribution comparison method. *Med Phys* 2003; 30: 2455–2464.
26. 日本放射線腫瘍学会 QA 委員会. 強度変調放射線治療における物理・技術的ガイドライン 2011
27. Alber M, Broggi S, De Wagter C, et al. Guidelines for the verification of IMRT, ESTRO Booklet No. 9. ESTRO, 2008, 89–106.
28. Cunningham JR, Cohen M, Walstam R, et al. Errors in dosimetry. In: ICRU REPORT 24, Determination of Absorbed Dose in a Patient Irradiated by Beams of X or Gamma Rays in Radiotherapy Procedures. ICRU, 1976, 45–50.

29. Van Dyk J and Purdy JA. Clinical implementation of technology and the quality assurance process. In: Van Dyk J eds. The modern technology of radiation oncology. Wisconsin, Medical Physics Publishing, 1999, 19-51.
30. Nelson CL, Mason BE, Robinson RCJ, et al. Commissioning results of an automated treatment planning verification system. *J Appl Clin Med Phys* 2014; 15: 4838.
31. Sun Nuclear Corporation. DoseCHECKTM & PerFRACTIONTM -On the Accuracy of the SNC Dose Calculator Algorithm. [https://www.sunnuclear.com/uploads/documents/whitepapers/DC\\_PF\\_Accuracy-of-the-SNC-Dose-Calculator\\_091216.pdf](https://www.sunnuclear.com/uploads/documents/whitepapers/DC_PF_Accuracy-of-the-SNC-Dose-Calculator_091216.pdf), (2 Sep 2022, date last accessed).
32. Narayanasamy G, Zalman T, Ha CS, et al. Evaluation of Dosimetry Check software for IMRT patient-specific quality assurance. *J Appl Clin Med Phys* 2015; 16: 5427.
33. Hoffmann L, Alber M, Söhn M, et al. Validation of the Acuros XB dose calculation algorithm versus Monte Carlo for clinical treatment plans. *Med Phys* 2018; 45: 3909-3915.
34. Kodama T, Saito Y, Hatanaka S, et al. Commissioning of the Mobius3D independent dose verification system for TomoTherapy. *J Appl Clin Med Phys*. 2019; 20: 12-20.
35. Jolly D, Dunn L, Kenny J. A clinical database to assess action levels and tolerances for the ongoing use of Mobius3D. *J Appl Clin Med Phys* 2017; 18: 59-65.
36. Ahmed S, Hunt D, Kapatoes J, et al. Validation of a GPU-Based 3D dose calculator for modulated beams. *J Appl Clin Med Phys* 2017; 18: 73-82.
37. Piffer S, Casati M, Marrazzo L, et al. Validation of a secondary dose check tool against Monte Carlo and analytical clinical dose calculation algorithms in VMAT. *J Appl Clin Med Phys* 2021; 22: 52-62.
38. Shepard AJ, Frigo SP. Characteristics and limitations of a secondary dose check software for VMAT plan calculation. *J Appl Clin Med Phys* 2021; 22: 216-223.
39. Stasinou D, Patatoukas G, Kollaros N, et al. Implementation of TG-218 for patient-specific quality assurance tolerance and action limits determination: Gamma passing rate evaluation using 3DVH software. *Med Phys* 2022; 49: 4322-4334.
40. 日本磁気共鳴医学会, 日本医学物理学会, 日本医学放射線学会, 日本放射線技術学会, 日本放射線腫瘍学会: MR 画像誘導即時適応放射線治療ガイドライン 2021 年版.
41. Zygmanski P, Rosca F, Kadam D, et al. Determination of depth and field size dependence of multileaf collimator transmission in intensity-modulated radiation therapy beams. *J Appl Clin Med Phys* 2007; 8: 76-95.
42. Vieilleveigne L, Khamphan C, Saez J, et al. On the need for tuning the dosimetric leaf gap for stereotactic treatment plans in the Eclipse treatment planning system. *J Appl Clin Med Phys* 2019; 20: 68-77.
43. Hernandez V, Hansen CR, Widesott L, et al. What is plan quality in radiotherapy? The importance of evaluating dose metrics, complexity, and robustness of treatment plans. *Radiother Oncol* 2020; 153: 26-33.

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1. 日本放射線腫瘍学会. 強度変調放射線治療における物理・技術的ガイドライン 2011.
2. Ding GX, Alaei P, Curran B et al. Image guidance doses delivered during radiotherapy: Quantification, management, and reduction: Report of the AAPM Therapy Physics Committee Task Group 180. *Med Phys*. May;45 2018.
3. 放射線医療技術書叢書 (33. 「放射線治療における位置照合とセットアップの実際」)
4. Verellen D: Image Guided Patient Setup. In: Bortfeld T eds. Image-guided IMRT. Springer, 97-116 2005.
5. Stanley DN, McConnell KA, Kirby N et al. Comparison of initial patient setup accuracy between surface imaging and three-point localization: A retrospective analysis. *J Appl Clin Med Phys* 2017; 18: 58-61.

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